



LiNA SeaStar™

REF: 4000-S-SL-NS, 4000-S-VL-NS, 4000-S-VS-NS,
4000-S-SS-NS, 4000-S-VV-NS, 4000-S-Q-NS

Caution: Federal Law (USA) restricts this device to sale by or on the order of a trained physician.

LiNA SeaStar™

REF: 4000-S-SL-NS, 4000-S-VL-NS, 4000-S-VS-NS,
4000-S-SS-NS, 4000-S-VV-NS, 4000-S-Q-NS

GB

Intended Use:

The LiNA SeaStar system, consisting of the SeaStar retractor ring and skin hooks, is intended to retract soft tissue from the operative surface.

Patient Population:

All groups dedicated to surgical procedures which require optimal exposure.

Precautions:

- For single use only. Do NOT reuse, reprocess or re-sterilize this device. Any reprocessing may impede the functions of this device. Reusing single use devices may also increase the risk of cross contamination. Attempts to clean the device results in risk of device malfunction.
- Use of the LiNA SeaStar requires adequate training and experience in performing surgical procedures. Intended user profile is a trained surgeon or physician.

Description of use:

Place the LiNA SeaStar around the initial incision and insert the hooks in directly opposing quadrants with equal tension. As dissection progresses, move the position of the hook into the deeper layers of the fascia. You may then add new hooks or reposition those already in place to achieve full balanced retraction from any angle. To ease closing maintain constant opposing tension.

Reporting:

Any serious incident that has occurred in relation to the device should be reported to the LiNA Medical ApS and the competent regulatory authority of the country in which the user and/or patient is established.

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