
# LiNA Bipolar Loop™

**Intended Use:**

The LiNA Bipolar Loop™ is a 5mm bipolar electrosurgical device intended for amputating the mobilized uterus during laparoscopic supracervical (subtotal) hysterectomy and resection of devascularized subserosal pedunculated myomas. To be used with an electrosurgical generator.

**Potential complications:**

- conversion to laparotomy;
- electrosurgical complications (e.g., insulation failure, electrocution);
- hemorrhage;
- implantable electrical device malfunction;
- small or large bowel injury, including perforation;
- urinary tract injury;
- thermal burns.

**Patient Population:**

Female of any age and weight dedicated to laparoscopic supracervical hysterectomy or subserosal pedunculated myomectomy.

**Contraindications:**

- The LiNA Bipolar Loop™ should **NOT** be used in the following situations:
- If the uterine isthmus cannot be visualized circumferentially before activation.
  - If the uterus or myoma is adhered to any adjacent organs.
  - If the patient has a large posterior uterine or cervical wall myomas.
  - Amputation of uterus which is larger than a gestational week of 20.

**Warnings:**

- The LiNA Bipolar Loop™ is supplied STERILE using an ethylene oxide process. Carefully inspect the packaging for any damage prior to use. Do NOT use if sterile barrier is damaged. Do NOT use past expiration date.
- For single use ONLY. Do NOT reuse, reprocess or re-sterilize the LiNA Bipolar Loop™. Any reprocessing may impede the functions of the device. Reusing single use devices may also increase the risk of cross contamination. Attempts to clean the device results in risk of device malfunction and/or erroneous pathology specimen collection due to residual tissue in the LiNA Bipolar Loop™.
- Do NOT use the device when in the presence of flammable gases or flammable fluids.
- Do NOT use if the LiNA Bipolar Loop™ is bent, twisted, malformed, or does not unfold completely to the original rhomboid design when the handle is completely advanced.
- Do NOT modify the original rhomboid design and electrode insulation of the LiNA Bipolar Loop™.
- The LiNA Bipolar Loop™ must be large enough to pass around the corpus of the uterus or the myoma. Then apply the loop around the cervix or the myoma stalk.
- Always verify that the active (i.e. un-insulated) part of the electrode is not in contact with any organ other than the uterus or myoma (e.g. bowel or vessel) before activating the generator.
- Do NOT activate the LiNA Bipolar Loop™ if it is not possible to visualize the cutting loop.
- **REMOVE the uterus manipulator before activating the LiNA Bipolar Loop™.** Only a uterus manipulator made of NON-CONDUCTING material (e.g. plastic) should be used with the LiNA Bipolar Loop™.
- **REMOVE intra uterine device (IUD) from the uterus prior to activating the LiNA Bipolar Loop™.**
- The LiNA Bipolar Loop™ should NOT be used in patients where the bladder cannot be completely dissected at the lower anterior uterine segment due to the risk of bladder injury.
- Do NOT allow the cable of the LiNA Bipolar Loop™ to be in contact with skin of the patient or surgeon during electrosurgical activation.
- Do NOT use excessive force/load to the LiNA Bipolar Loop™. This could damage the product or result in strangulation and mechanical sectioning of the uterus or myoma with the risk of excessive bleeding.
- Be sure to activate the electrode prior to retracting the LiNA Bipolar Loop™. Failure of activation can result in mechanical sectioning of the uterus or myoma with the risk of excessive bleeding.
- Do not use the LiNA Bipolar Loop™ for slicing the uterus or myoma as alternative to morcellation.
- Be sure that the ceramic tip is fully retracted into the shaft /tube prior to extracting the LiNA Bipolar Loop™ through the trocar to prevent the ceramic tip from getting caught in the flap of the valve.
- After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.
- DO NOT activate if the uninsulated wires of LiNA Bipolar Loop™ are in contact with metal objects.

**Precautions:**

- Use of the LiNA Bipolar Loop™ requires adequate training and experience in performing laparoscopic hysterectomy and myomectomy.
- Failure to carefully follow all applicable instructions may result in significant injury to the patient, physician or attendants and may have an adverse effect on the outcome of procedures performed.

- Precautions should be taken when the patient has any implanted electronic device (IED). The perioperative team should consult with the team managing the implanted device before surgery and as needed. When LiNA Bipolar Loop™ is activated it should be kept more than ½ inch (1 cm) from the IED.
- Carefully use the LiNA Bipolar Loop™ with other sharp tools as wire's insulation damage can lead to device malfunction.

**Product description:**

The LiNA Bipolar Loop™ consists of a black outer-tube and a blue inner-tube at which the handle of the device is attached. The rhomboid loop is fastened to the inner-tube and is advanced by holding the upper part of the black outer-tube while pushing forward the blue handle. To retract the loop, gently pull the blue handle with one hand while holding the black outer-tube with the other.

**Generator specifications & settings:**

The LiNA Bipolar Loop™ is compatible with most standard electrosurgical generators with bipolar outlet (fixed 2-pin 4mm plug with a 29 mm distance).

- The recommended setting for the generator is 55-70 Watts. Use the lowest possible settings capable of achieving the desired surgical effect.
- Use a bipolar foot pedal.
- Use only electrosurgical generators working in the frequency field 350 to 472 KHz as a minimum covering the power range 55-70 Watts of bipolar energy.

**Instructions for use:**

This information is provided for using the LiNA Bipolar Loop™. It is not intended to demonstrate laparoscopic techniques or procedures. Consult published and recognized medical literature for appropriate techniques, complications, and hazards prior to doing any laparoscopic procedure. Consult the instructions for use provided by the manufacturer of the electrosurgical generator and/or laparoscope to be used for contraindications. Read the instruction for use prior to use of this device!

**After the uterine ligaments and vessels have been correctly ligated, the uterus can be amputated from the cervix or the myoma from the stalk, following the instructions below:**

1. Carefully inspect the packaging for any damage prior to use. Do NOT attempt to use the device if sterile barrier is damaged.
2. Connect the LiNA Bipolar Loop™ to the generator's bipolar outlet. Set the generator between 55-70 Watts depending on type of generator.
3. Carefully introduce the LiNA Bipolar Loop™ through a 5mm trocar.
4. **a.** Unfold the electrode (loop) by forwarding the handle. The handle must be completely advanced and the loop 100% unfolded prior to placing the loop around the uterine body or myoma. **b.** Use grasper(s) and the uterine manipulator to assist placing the electrode (loop) around the uterine body or myoma.
5. Before activating the generator always verify that the active (i.e. un-insulated) portion of the electrode is not in contact with any organ (e.g. bowel or vessel) other than the uterus or myoma.
6. **a.** Carefully tighten the loop around the lower part of uterine body (above cervix) or myoma stalk. **b. REMOVE** the uterine manipulator by slightly releasing the tension of the loop.
7. **a.** Activate the LiNA Bipolar Loop™ by pressing the bipolar pedal. Be sure to await the bipolar energy before sectioning the uterus or myoma. Immediate retraction of the loop may inflict an unintended mechanical cut. **b.** After amputating the uterus or myoma, do not activate the electrodes to manipulate or divide the uterus or myoma.
8. After amputating the uterus from the cervix, ensure that the proximal portion of the endocervical canal which may contain endometrial cells is coagulated, and ensure that any bleeding or oozing on the cervical stump is coagulated.
9. Remove the LiNA Bipolar Loop™ through the trocar.
10. Dispose the LiNA Bipolar Loop™ after usage.

Provide instructions for follow-up to patients in the event of complications including intra abdominal bleeding which in rare cases can lead to exsanguinations.


**Reporting:**

Any serious incident that has occurred in relation to the device should be reported to the LiNA Medical ApS and the competent regulatory authority of the country in which the user and/or patient is established.

**Caution:**

Federal Law (USA) restricts this device to sale by or on the order of a trained physician.

Country of origin: Poland.  
CE marking 2012.

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