LiNA Skin Hook™ (sterile)

Ref: SH-110, SH-120, SH-220, SH-130

en

Intended use:

The LiNA Skin Hook™ is intended to retract the soft tissue from the operative surface.

Patient population:

Patients intended to have surgical procedures which require optimal exposure.

Product description:

The LiNA Skin Hook™ together with LiNA SeaStar™ is a complete, self-retaining system for retracting delicate soft tissue. By securing the LiNA Skin Hook™ on to the LiNA SeaStar™ the exposure of the incision is increased. Den fleksible og elastiske haakke kan justeres under the procedure to create optimum tension.

Precautions:

- For single use only. Do not reuse, reprocess or re-sterilize this device. Any reprocessing may impede the functions of this device. Reusing single use devices may also increase the risk of cross contamination. Attempts to clean the device results in risk of device malfunction.
- Use of the LiNA Skin Hook™ requires adequate training and experience in performing surgical procedures. Intended user profile is a trained surgeon or physician.
- After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy. Be aware of the sharp edges.
- Please do not attach LiNA Skin Hook™ to home-made retractor rings.

Warnings:

- The LiNA Skin Hook™ is supplied sterile using an ethylene oxide process. Carefully inspect the packaging for any damage prior to use. Do not use if the sterile barrier is damaged. Do not use past expiration date.
- Carefully inspect the device for any damage prior to use. Dispose of the device if bent, malformed or if breaks or cracks are found.
- Special attention is needed when installing and removing LiNA Skin Hook™ from the operative surface since device has sharp edges.

Instruction for use:

- Remove LiNA Skin Hook™ from the pouch. Gently slide the silicone tube from the hook starting from the end of the tube. Do not use excessive force as it may cause product damage or user injury.
- Place the LiNA SeaStar™ around the initial incision and insert the LiNA Skin Hook™ in directly opposing quadrants with equal tension.
- As dissection progresses, move the position of the hook into the deeper layers of the fascia.
- You may then add new hooks or reposition those already in place to achieve full balanced retraction from any angle.
- To ease closing maintain constant opposing tension.

Reporting:

Any serious incident that has occurred in relation to the device should be reported to the LiNA Medical ApS and the competent regulatory authority of the country in which the user and/or patient is established.

Caution:

Federal Law (USA) restricts this device to sale by or on the order of a trained physician.

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