LiNA Gold Loop™

STERILE EO

CE
0459











LiNA Gold Loop™

REF EL-160-4, EL-160-8, EL-200-4, EL-200-8, EL-240-4, EL-240-8

LiNA Gold Loop™

REF EL-160-4, EL-160-8, EL-200-4, EL-200-8, EL-240-4, EL-240-8

Important:

This information is provided for using the LiNA Gold Loop™. It is not intended to demonstrate laparoscopic techniques or procedures. Consult published and recognized medical literature for appropriate techniques, complications, and hazards prior to doing any laparoscopic procedure. Consult the instructions for use provided by the manufacturer of the electrosurgical generator and/or laparoscope to be used for contraindications.

Intended use:

The LiNA Gold Loop™ is a monopolar device intended for Laparoscopic Supracervical Hysterectomy (LSH) and Subserosal Pedunculated Myomectomy.

Patient population:

Intended patient population are women undergoing Laparoscopic Supracervical Hysterectomy (LSH) and Subserosal Pedunculated Myomectomy.

Product description:

The LiNA Gold Loop™ is a 5mm monopolar electrosurgical device. It consists of wire loop inside an introducer sheath which can be inserted into the abdominal cavity through a laparoscopic access port. The loop opens when the handle is pushed forward. Rhomboid shaped loop is available in different sizes depending on LiNA Gold Loop™ variant.

Contraindications:

The LiNA Gold Loop™ should not be used in the following situations:

- If the uterine isthmus cannot be visualized circumferentially before activation.
- If the uterus is adhered to any adjacent organs.
- If the patient has a large posterior uterine or cervical wall myomas.
- Removal of uterus which are larger than a gestational week of 20.
- Removal of uterus which is attached to the pelvis with large laterally fibroids.
- The LiNA Gold Loop™ should not be used in patients who have undergone caesarean delivery due to the risk of bladder injury in the event that the bladder cannot be completely dissected off of the lower anterior uterine segment.

Warnings:

- The LiNA Gold Loop™ is supplied sterile using an ethylene oxide process. Do not use if sterile barrier is damaged.
- Carefully inspect the packaging for any damage prior to use. Do not use if the packaging is damaged.
- For single use only. Do not reuse, reprocess or re-sterilize this device. Any reprocessing may impede the functions of this device. Reusing single use devices may also increase the risk of

cross contamination. Attempts to clean the device results in risk of device malfunction.

- Carefully inspect the device for unintended rough surfaces, sharp edges, protrusions or other damages prior to usage. Make sure the shaft is not bended. If any damage is discovered, dispose of the device.
- Do not use the LiNA Gold Loop™ when in the presence of flammable material or flammable gasses/fluids.
- Do not allow the cable of the LiNA Gold Loop™ to be in contact with skin of the patient or surgeon during electrosurgical activation.
- The LiNA Gold Loop™ shall be used with monopolar electrosurgical units, with a maximum operating voltage not exceeding 2.5kV peak, and may only be used with a EN-60601-1 approved generator.
- The LiNA Gold Loop™ must be large enough to pass around the corpus of the uterus. Then apply the loop around the cervix.
- The LiNA Gold Loop™ must not be pre-tested up against bodily tissue prior to any cutting process.
- Always verify that the active (i.e. un-insulated) portion of the electrode is not in contact with any organ other than the uterus (e.g. bowel or vessel) before activating the generator.
- Do not use excessive force while using LiNA Gold Loop™. This could damage the product or result in entrapment and mechanical sectioning of the uterus or myoma with the risk of excessive bleeding.
- Do not use the LiNA Gold Loop™ for slicing the uterus or myoma as alternative to morcellation.
- Do not activate the LiNA Gold Loop™ if it is not possible to visualize the cutting loop.
- Remove the uterus manipulator before activating the LiNA Gold Loop™. Only uterus manipulator made of non conducting material (e.g. plastic) should be used with the LiNA Gold Loop™.
- Do not use if the loop is bent, twisted, malformed, or does not unfold completely to the original rhomboid design when the handle is completely advanced.
- Do not modify the original rhomboid design of the LiNA Gold Loop™.
- No modification of this equipment is allowed.
- After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

Precautions:

- Use of the LiNA Gold Loop™ requires adequate training and experience in performing laparoscopic myomectomy and hysterectomy.
- The user must be trained and competent in the basics of HF surgery and safety precautions that apply to the use of HF surgery.
- Precautions should be taken when the patient has any implanted electronic device (IED).

Generator specifications and settings:

- The LiNA Gold Loop™ is compatible with all standard electrosurgical generators with monopolar outlet and BLEND modes.
- The recommended setting for the generator is 100 - 130 Watt, Blend 1. Max. voltage 2500 Vp-p.

Important:

- Keep the voltage/power as low as possible to achieve the desired effect.
- Grounding plate is necessary.
- This device should only be used by persons qualified in the procedures for which this is indicated.

Instruction for use:

After the uterine ligaments and vessels have been correctly ligated, when the uterus can be safely amputated from the cervix, follow the instructions below:

1. Apply the grounding plate to the patient and generator. Connect the LiNA Gold Loop™ to the generator's monopolar outlet. Set the generator at Blend 1 and start at 100 Watt.
2. Carefully insert the LiNA Gold Loop™ through a 5mm trocar/entry.
3. Unfold the electrode (loop) by forwarding the handle. The handle must be completely advanced and the loop 100% unfolded prior to placing the loop around the uterine body.
- 4a. Use grasper(s) and the uterine manipulator to help place the electrode (loop) around the uterine body.

4b. Before activating the generator always verify that the active (i.e. un-insulated) portion of the electrode is not in contact with any organ other than the uterus (e.g. bowel or vessel).

5a. Carefully tighten the loop around lower part of uterine body (above cervix).

5b. Remove the uterine manipulator and slightly release the tension of the loop.

6. Start sectioning the uterine body by pressing "cut" on the foot pedal. After 1-3 seconds slowly pull the handle, until the electrode (loop) has cut the uterine body and is retracted inside the tube. Stop pressing the foot pedal.

7. Remove the LiNA Gold Loop™ through the trocar.

8. Dispose of the used LiNA Gold Loop™ after usage.

- After amputating the uterus from the cervix, ensure that the proximal portion of the endocervical canal which may contain endometrial cells is coagulated, and ensure that any bleeding or oozing on the cervical stump is coagulated.
- Provide instructions for follow up to patients in the event of complications including intra abdominal bleeding which in rare cases can lead to exsanguinations.

Reporting:

Any serious incident that has occurred in relation to the device should be reported to the LiNA Medical ApS and the competent regulatory authority of the country in which the user and/or patient is established.

Country of origin: Poland.

CE marking 2006.

Caution: Federal Law (USA) restricts this device to sale by or on the order of a trained physician.

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