
# LiNA SeaStar™

(Sterile)

**Intended use:**

The LiNA SeaStar™ is intended to retract soft tissue from the operative surface.

**Patient population:**

All groups dedicated to surgical procedures which require optimal exposure.

**Device description:**

The LiNA SeaStar™ is a retractor ring with grooves every 2,5 cm on the convex site. The hinges are controlled by finger screws allowing the two ring halves to remain at any desired angle. The device is provided in a broad range of shapes and sizes which can easily be adjusted to fit various anatomical areas, thereby facilitating different surgical applications. The LiNA SeaStar™, together with LiNA Skin Hook™ is a complete system for retracting delicate soft tissue. The flexible design of the system allows easy adjustment and repositioning of the LiNA SeaStar™ during the surgical procedures and provides optimal exposure of the surgical site.

**Precautions:**

- For single use only. Do NOT reuse, reprocess or re-sterilize this device. Any reprocessing may impede the functions of this device. Reusing single use devices may also increase the risk of cross contamination. Attempts to clean the device results in risk of device malfunction.
- The LiNA SeaStar™ is supplied sterile using an ethylene oxide process. Carefully inspect the packaging for any damage prior to use. Do not use if sterile barrier is breached.
- No modification of this equipment is allowed.
- Do not use excessive force during the screw tightening.
- Use of the LiNA SeaStar™ requires adequate training and experience in performing surgical procedures. Intended user profile is a trained surgeon or physician.
- After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

**Warnings:**

- Do not use past expiration date.
- Carefully inspect the device for any damage prior to use. Dispose of the device if bent, malformed or if breaks or cracks are found.

**Instruction for Use:**

Place the LiNA SeaStar™ around the initial incision and insert the hooks in directly opposing quadrants with equal tension. As dissection progresses, move the position of the hook into the deeper layers of the fascia. You may then add new hooks or reposition those already in place to achieve full balanced retraction from any angle. To ease closing maintain constant opposing tension.

**Reporting:**





Any serious incident that has occurred in relation to the device should be reported to the LiNA Medical ApS and the competent regulatory authority of the country in which the user and/or patient is established.

**Caution:**

Federal Law (USA) restricts this device to sale by or on the order of a trained physician.

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