










LiNA Bipolar Loop™

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

en

LiNA Bipolar Loop™

REF BL-160-E, BL-200-E, BL-240-E

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

Product description:

The LiNA Bipolar Loop™ consists of a black outer-tube and a blue inner-tube at which the handle of the device is attached. The rhomboid loop is fastened to the inner-tube and is advanced by holding the upper part of the black outer-tube while pushing forward the blue handle. To retract the loop, gently pull the blue handle with one hand while holding the black outer-tube with the other.

Intended Use:

The LiNA Bipolar Loop™ is a 5mm bipolar electro-surgical device intended for amputating the mobilized uterus during laparoscopic supracervical (subtotal) hysterectomy and resection of devascularized subserosal pedunculated myomas.

Indication for use:

The LiNA Bipolar Loop™ is a 5mm bipolar electro-surgical device intended for amputating the mobilized uterus during laparoscopic supracervical (subtotal) hysterectomy and resection of devascularized subserosal pedunculated myomas. To be used with an electro-surgical generator that provides a bipolar outlet.

Patient population:

Women requiring laparoscopic supracervical hysterectomy or subserosal pedunculated myomectomy.

Contraindications:

The LiNA Bipolar Loop™ should not be used in the following situations:

- If the uterine isthmus cannot be visualized circumferentially before activation.
- If the uterus or myoma is adhered to any adjacent organs.
- If the patient has large posterior uterine or cervical wall myomas.
- Amputation of uterus which is larger than a gestational week of 20.

Potential complications:

- Unintended burns (might lead to hemorrhage, small or large bowel injury including perforation, urinary tract injury, conversion to laparotomy),
- Implantable electrical device malfunction,
- Electrical shock,
- Complications related with surgical smoke,
- Explosion,
- Surgical fire,
- Electromagnetic interference with other devices,
- Bleeding,
- Infection,
- Pain.

Warnings:

- The LiNA Bipolar Loop™ is provided sterile using ethylene oxide sterilization. Carefully inspect the packaging for any damage prior to use. Do not attempt to use the device if product sterile barrier or its packaging is damaged.
- For single use only. Do not reuse, reprocess or re-sterilize the LiNA Bipolar Loop™. Any reprocessing

may impede the functions of the device. Reusing single use devices may also increase the risk of cross contamination. Attempts to clean the device results in risk of device malfunction.

- Carefully inspect the device before use. If any damage or defect is identified, dispose of the device.
- No modification of this equipment is allowed. Do not attempt to modify the original rhomboid design and electrode insulation of the LiNA Bipolar Loop™.
- Do not use if the LiNA Bipolar Loop™ is bent, twisted, malformed, or does not unfold completely to the original rhomboid design when the handle is completely advanced.
- Do not use past expiration date.
- Do not use if the device is exposed to non-sterile surfaces before procedure.
- Do not use in the presence of flammable anesthetics or oxidizing gases (such as nitrous oxide (N₂O) and oxygen) or in close proximity to volatile solvents (such as ether or alcohol), as explosion may occur.
- Do not place instruments near or in contact with flammable materials (such as gauze or surgical drapes). Instruments that are activated or hot from use may cause a fire.
- The LiNA Bipolar Loop™ must be large enough to pass around the corpus of the uterus or the myoma. Then apply the loop around the cervix or the myoma stalk.
- The LiNA Bipolar Loop™ shall be used with bipolar electro-surgical units, with a maximum operating voltage not exceeding 1.4kV peak, and may only be used with a IEC 60601-1, IEC 60601-2-2 and IEC 60601-1-2 approved generator.
- Always verify that the active (i.e. un-insulated) part of the electrode is not in contact with any organ other than the uterus or myoma (e.g. bowel or vessel) before activating the generator.
- Avoid direct contact between the active electrode and other metal surgical equipment during surgery as this can result in thermal damage to tissue. Always keep the active electrode in view during the procedure.
- The LiNA Bipolar Loop™ should not be used in patients where the bladder cannot be completely dissected at the lower anterior uterine segment due to the risk of bladder injury.
- Do not allow the cable of the LiNA Bipolar Loop™ to be in contact with skin of the patient or surgeon during electro-surgical activation.
- Do not use excessive force/load to the LiNA Bipolar Loop™. This could damage the product or result in strangulation and mechanical sectioning of the uterus or myoma with the risk of excessive bleeding.
- Remove the uterus manipulator before activating the LiNA Bipolar Loop™. Only a uterus manipulator made of non-conducting material (e.g. plastic) should be used with the LiNA Bipolar Loop™.
- Remove intra uterine device (IUD) from the uterus prior to activating the LiNA Bipolar Loop™.
- Be sure to activate the electrode prior to retracting the LiNA Bipolar Loop™. Failure of activation can result in mechanical sectioning of the uterus or myoma with the risk of excessive bleeding.
- Do not use the LiNA Bipolar Loop™ for slicing the uterus or myoma as alternative to morcellation.
- Be sure that the ceramic tip is fully retracted into the shaft /tube prior to extracting the LiNA Bipolar Loop™

LiNA Bipolar Loop™

REF BL-160-E, BL-200-E, BL-240-E

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

through the trocar to prevent the ceramic tip from getting caught in the flap of the valve.

- When not using LiNA Bipolar Loop™, place it in a clean, dry, highly visible area not in contact with the patient. Inadvertent contact with the patient may result in burns.
- Due to concerns about the carcinogenic and infectious potential of electro-surgical byproducts (such as tissue smoke plume and aerosols), protective eyewear, filtration masks, and effective smoke evacuation equipment should be used.
- Pay attention when the patient has any implantable electronic device (IED). The perioperative team shall consult with the team managing the implanted device before surgery and as needed. In case team managing the device is not available, consult with the device manufacturer or interdisciplinary team in your hospital. When LiNA Bipolar Loop™ is activated it shall be kept more than ½ inch (1 cm) from the IED.
- Failure to follow all instructions or any warnings or precautions may result in significant injury to the patient, physician or attendants and may have an adverse effect on the outcome of procedures performed. LiNA Bipolar Loop™ should not be used for any other purpose than intended function.
- Connect the LiNA Bipolar Loop™ to the electro-surgical unit only when the energy is off. Failure to do so may result in an injury or electrical shock to the patient or operating room personnel.

Precautions:

- Always have a backup device readily available for immediate use. If any malfunction should occur during use, stop the procedure immediately, and slowly withdraw LiNA Bipolar Loop™ and replace with a new device.
- Use of the LiNA Bipolar Loop™ requires adequate training and experience in performing laparoscopic hysterectomy and myomectomy.
- Carefully use the LiNA Bipolar Loop™ with other sharp tools as wire's insulation damage can lead to device malfunction.
- After use, the device is to be considered contaminated medical waste and may pose a potential infection risk. Dispose of the device and packaging in accordance with applicable laws, administrative, and/or local government policy and hospital procedures, including those regarding biohazards, microbiological hazards and infectious substances.

Generator specifications and settings:

The LiNA Bipolar Loop™ is compatible with most standard electro-surgical generators with bipolar outlet. Connect the European connector to a standard bipolar cable which fits your generator. Connect the bipolar cable to the bipolar outlet on the generator.

- The recommended setting for the generator is 55-70 Watts. Use the lowest possible settings capable of achieving the desired surgical effect.
- Use a bipolar foot pedal.
- Use only electro-surgical generators working in the frequency field 350 to 472 KHz as a minimum covering the power range 55-70 Watts of bipolar energy.

Instruction for use:

After the uterine ligaments and vessels have been correctly ligated, the uterus can be amputated from the cervix or the myoma from the stalk, following the instructions below:

1. Connect the LiNA Bipolar Loop™ to the generator's bipolar outlet. Set the generator between 55-70 Watts depending on type of generator.
2. Carefully introduce the LiNA Bipolar Loop™ through a 5 mm trocar.
3. a. Unfold the electrode (loop) by forwarding the handle. The handle must be completely advanced and the loop 100% unfolded prior to placing the loop around the uterine body or myoma.
b. Use grasper(s) and the uterine manipulator to assist placing the electrode (loop) around the uterine body or myoma.
4. Before activating the generator always verify that the active (i.e. un – insulated) portion of the electrode is not in contact with any organ (e.g. bowel or vessel) other than the uterus or myoma.
5. a. Carefully tighten the loop around the lower part of uterine body (above cervix) or myoma stalk.
b. REMOVE the uterine manipulator by slightly releasing the tension of the loop.
6. a. Activate the LiNA Bipolar Loop™ by pressing the bipolar pedal. Be sure to await the bipolar energy before sectioning the uterus or myoma. Immediate retraction of the loop may inflict an unintended mechanical cut.
b. After amputating the uterus or myoma, do not activate the electrodes to manipulate or divide the uterus or myoma.
7. After amputating the uterus from the cervix, ensure that the proximal portion of the endocervical canal which may contain endometrial cells is coagulated, and ensure that any bleeding or oozing on the cervical stump is coagulated.
8. Remove the LiNA Bipolar Loop™ through the trocar.
9. Dispose the LiNA Bipolar Loop™ after usage. Provide instructions for follow-up to patients in the event of complications including intra-abdominal bleeding which in rare cases can lead to exsanguinations.

Reporting:

Any serious incident that has occurred in relation to the device should be reported to the LiNA Medical AG and the competent regulatory authority of the country in which the user and/or patient is established.

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