

Cystoscopy without  
capital cost or  
complexity



## LiNA OperåScope™

Single-use Cystoscopy System

The complete, portable system  
for cystoscopy.



# LiNA OperåScope™

Single-use Cystoscopy System

The LiNA OperåScope™ single-use cystoscopy system is the first and only fully disposable operative cystoscope. Designed for use in visualization of the female urethra and the bladder cavity during diagnostic and therapeutic cystoscopic procedures.

- **Convenient solution** without capital cost or complexity
- **Single-use** design is ready for use out of the box
- **360° steerable, deflected tip** allows for precise control and maneuverability
- **Slim 4.2mm catheter** designed for improved patient tolerability
- **5.5 french operative channel and single-use instrumentation** for therapeutic use

## Ordering information:

Product Code	Description	UOM
OP-201-6	LiNA OperåScope™ Operative Hysteroscope with HDMI cable and onboard LCD (sterile single-use)	Box of 6
OP-FCP-6	LiNA OperåScope™ Biopsy Forceps (sterile single-use)	Box of 6
OP-GRA-6	LiNA OperåScope™ Rat Tooth Alligator Grasper (sterile single-use)	Box of 6
OP-RM-1	LiNA OperåScope™ Recording Module	Each
OP-TDK-6	LiNA OperåScope™ Tubing & Drape Kit (sterile single-use)	Box of 6

### LiNA Medical ApS

Formervangen 5  
DK-2600 Glostrup, Denmark  
Tel.: +45 43 29 66 66

Fax: +45 43 29 66 99  
info@lina-medical.com  
www.lina-medical.com  
CVR: 10054974

### Distributor:

**LiNA Medical USA**  
Norcross, GA 30093, USA

Tel: +1 (855) 546-2633  
Fax: +1 (678) 380-0500  
info@linamed.com  
www.linamed.com

©2019 LiNA Medical ApS. All rights reserved. All products may not be available for sale in all geographies. Please contact your local LiNA Medical representative for information about products available in your area. This information is intended for healthcare professionals. If you are a patient seeking more information, please consult your healthcare provider. Before using any LiNA product, always refer to the IFU included in our packaging for complete indications for use, contraindications, warnings, precautions, and instructions. LiNA is a registered trademark of LiNA Medical in EU, US and CN.

