

Cystoscopy without
capital cost or
complexity



LiNA Operascope™

Single-use portable cystoscope



LiNA OperåScope™

Single-use portable cystoscope


The LiNA OperåScope™ single-use portable cystoscope is the first and only fully disposable operative cystoscope. Designed for use in visualization of the female urethra and the bladder cavity during diagnostic and therapeutic cystoscopic procedures.

- **Convenient solution** without capital cost or complexity
- **Single-use** design is ready for use out of the box
- **360° steerable, deflected tip** allows for precise control and maneuverability
- **Slim 4.2mm catheter** designed for improved patient tolerability
- **5.5 french operative channel** for therapeutic use

Ordering information:

Product Code	Description	UOM
OP-201-6-I	LiNA OperåScope™ Operative Hysteroscope with HDMI cable and onboard LCD (sterile, single-use)	Box of 6
OP-RM-1-UK	LiNA OperåScope™ Recording Module – Socket Type G	Each
OP-RM-1-EU	LiNA OperåScope™ Recording Module – Socket Type C	Each
OP-RM-1-AUS	LiNA OperåScope™ Recording Module – Socket Type I	Each

LiNA OperåScope™



LiNA Medical ApS Formervangen 5, DK-2600 Glostrup, Denmark

Tel.: +45 43 29 66 66 | Fax: +45 43 29 66 99 | info@lina-medical.com | www.lina-medical.com | CVR: 10054974

©2020 LiNA Medical ApS. All rights reserved.
LiNA is a registered trademark of LiNA Medical in EU, US and CN.

EN ISO 13485:2016
Quality Management System

