


LiNA OperåScope® Needle

Single-use rigid needle, ideal for in-office cystoscopic injections



Innovation in Gynecology


LiNA OperåScope® Needle

Single-use rigid needle, ideal for in-office cystoscopic injections

Simplicity and Versatility

- Sterile, single-use
- Two sizes – to accommodate bladder wall thickness
- Protection cap
- Limit buckle to control needle exposure during introduction
- Luer lock for syringe connection

Available LiNA OperåScope® Needle sizes



Technical Specifications:

LiNA OperåScope® Needle	
Required channel inner diameter	5 Fr
Working length	313 mm
Needle width	23G
Needle length	4 or 8 mm
Tip shape	3-bevel tip

Ordering information:

Product code	Description	LiNA OperåScope™ Needle UOM
OP-201-6	LiNA OperåScope® with HDMI cable & LCD (sterile, single-use)	Box of 6 ea
OP-NEE4-6	LiNA OperåScope® Needle, 4 mm (sterile, single-use)	Box of 6 ea
OP-NEE8-6	LiNA OperåScope® Needle, 8 mm (sterile, single-use)	Box of 6 ea

LiNA Medical ApS
Formervangen 5
DK-2600 Glostrup, Denmark
Tel.: +45 43 29 66 66

Fax: +45 43 29 66 99
info@lina-medical.com
www.lina-medical.com
CVR: 10054974

Distributor:
LiNA Medical USA
Norcross, GA 30093, USA

Tel: +1 (855) 546-2633
Fax: +1 (678) 380-0500
info@linamed.com
www.linamed.com

©2023 LiNA Medical ApS. All rights reserved. All products may not be available for sale in all geographies. Please contact your local LiNA Medical representative for information about products available in your area. This information is intended for healthcare professionals. If you are a patient seeking more information, please consult your healthcare provider. Before using any LiNA product, always refer to the IFU included in our packaging for complete indications for use, contraindications, warnings, precautions, and instructions. LiNA is a registered trademark of LiNA Medical in EU, US and CN.

EN ISO 13485:2016
Quality Management System

