


LiNA OperåScope™ Lasso

For hysteroscopic polypectomies



Innovation in Gynecology

LiNA OperåScope™ Lasso

For hysteroscopic polypectomies

The LiNA OperåScope™ Lasso is a single use mechanical instrument, specifically designed for removing intrauterine polyps.

Versatile

- Available in two sizes


Ready to use, out of the box

- Sterile, single-use

No energy

- Mechanical instrument

Compatible with hysteroscopes with a straight working channel of minimum Ø1.87 mm diameter and a maximum length of 310 mm.



Ordering Information:

LiNA OperåScope™ and LiNA OperåScope™ Lasso



Product code	Description	UOM
OP-201-6-I	LiNA OperåScope™ with HDMI cable & LCD (sterile, single-use)	Box of 6ea
OP-LAS10-3-I	LiNA OperåScope™ Lasso 10mm Basket (sterile, single-use)	Box of 3ea
OP-LAS16-3-I	LiNA OperåScope™ Lasso 16mm Basket (sterile, single-use)	Box of 3ea

Ordering Information:

LiNA Versascope™



Product code	Description	UOM
GMS805	LiNA Versascope™ Sheath	Box of 5ea

Distributor



EU MDR 2017/745 device class for REF:
- OP-201-6-I class IIa
- OP-LAS10-3-I, OP-LAS16-3-I: class Is
Name and number of the Notified Body: TÜV Nord Polska Sp. z o.o. # 2274.
- GMS805: class Is
Name and number of the Notified Body: TÜV Rheinland LGA Products GmbH #0197

LiNA OperåScope Lasso / 2025-06-BV0269C

©2026 LiNA Medical AG. All rights reserved. All products may not be available for sale in all geographies. Please contact your local LiNA Medical representative for information about products available in your area. This information is intended for healthcare professionals. If you are a patient seeking more information, please consult your healthcare provider. Before using any LiNA product, always refer to the IFU included in our packaging for complete indications for use, contraindications, warnings, precautions, and instructions.

LiNA is a registered trademark of LiNA Medical in EU, US and CN.

LiNA Medical AG, Platz 3, 6039 Root D4 | Switzerland
Phone: +41 78 617 56 56, info@lina-medical.com | www.linamed.com

LiNA Medical Polska Sp. z o.o., Rolna 8A, Sady, 62-080 Tarnowo Podgórze | Poland
Phone: +48 61 222 21 21, info.poland@lina-medical.com | www.lina-medical.pl

